

800 SERIES


>>> Δυτοματις <

Call:+44 (0)1782 339320 | Email: sales@iae.co.uk | Weblink: https://iae.co.uk/




800 SERIES (B21) AUTO YOKE SQUEEZE CRUSH

INSTALLATION, OPERATION & MAINTENANCE TECHNICAL DOCUMENT NO. 1002240-01



PRODUCT 1002395-01 KEY FEATURES



INTENDED USE

A heavy-duty restraining device that holds big, heavy and strong animals like cattle so that work force (including veterinarians) can work with them. It is close fitting and simply keeps the animals from moving, which keeps both the cattle and workers safe.

TECHNICAL DATA

Item Number: 1002395-01 Load: 1000 Kg Max. Weight: 695 Kg Dimensions (W x H x D): 1225 x 2230 x 2365 Approx.

AUTO YOKE

Auto yoke is pre-set to 'ready to receive' mode by personnel and subsequently it is the animal that controls the yoking process by its forward movement - is designed for use in environments where cattle are frequently handled in high numbers

Auto yoke features include,

□ NEW DESIGN self-locking brake box mechanism

 \Box Reset bolt on operating handle \triangle also used for emergency yoke gate release \triangle

□ Yoke adjustment achievable, gap adjustment ranges from 240mm to 100mm (full height) spring-loaded yoke adjustment link arm

COMPATIBLE WITH IAE OPTIONS TO FIT

- CRUSH EXTENSION 1002145-01
- SLIDING GATE IN FRAME 1001686-01
- A800 SCOOP BRACKET 1003653-01
- 800 SERIES ID READER KIT 1003656-01
- HOOF TRIMMING KIT 1002393-01
- STERNUM BAR KIT 1002705-01
- LOAD BAR

(All options to fit supplied separate with onsite assembly required, to be organized by customer/operator)

SAFETY INFORMATION

PLEASE TAKE NOTE OF THE FOLLOWING SAFETY PRECAUTIONS BEFORE USE FOR THE WELLBEING OF ALL ANIMALS AND PERSONNEL. ANY HARM OR INJURY WHICH OCCURS FROM THE LACK OF COMPLIANCE WITH THESE HEALTH AND SAFETY MEASURES WILL BE THE RESPONSIBILITY OF THE OPERATOR.

- □ **IMPORTANT!** Read these instructions carefully before use and keep them for future reference. The safety of your animals may be affected if you do not follow these instructions in this manual.
- □ WARNING! DANGER TO LIFE AND RISK OF ACCIDENTS FOR INFANTS AND CHILDREN! Keep children away from the product when in operation.
- □ Situate this product on a level, solid surface.
- □ WARNING! Be careful with hands, fingers and limbs when operating the crush.
- Personal protection may be required for some stages of installation, and/or operation.
- □ **CAUTION!** When handling heavy equipment, some stages of installation may require at least two people.
- DO NOT lean on or exert undue pressure on the product.
- CAUTION! Check posture when lifting/manoeuvring the equipment or associated parts.
- □ ALWAYS Check the product for damage before each use.
- □ NEVER USE THE PRODUCT IF IT IS DAMAGED.

If you have any questions or comments regarding the product, we will be glad to hear them either through our website www.iae.co.uk or by contacting your local IAE dealership.

INTRODUCTION/PRODUCT FEATURES	Page 2	
INTENDED USE	Page 3	
TECHNICAL DATA	Page 3	
COMPATIBLE WITH IAE OPTIONS TO FIT	.Page 3	
SAFETY INFORMATION	Page 4	
CRUSH INSTALL/GROUND FIXING	Page 6	I
HOW TO OPERATE		
OPERATING AUTO YOKE	Page 7	
YOKE APERTURE ADJUSTMENT GUIDE	Page 9	
SQUEEZE MECHANISM	Page 1	0
OPERATING SIDE DOORS	Page 1	1
OPERATING FULL LENGTH SIDE RAIL	Page 1	1
OPERATING VET ACCESS PANELS	Page 1	2
REMOVING VET ACCESS PANELS	Page 1	2
PIVOTING VET ACCESS TOP PANEL	Page 1	3
OPERATING RUMP RAIL	.Page 1	4
RECOMMENDED MAINTENANCE	Page 1	6
CONTACT	Page 1	6


CRUSH FIRST USE/GROUND FIXING

- Locate on firm and level ground, IAE recommend placing onto a concrete pad.
- Allow for sufficient space all around the crush for safe operation.
- It is advised to run through the operations of the crush without cattle present to familiarise yourself with the correct workings.

Ground fixing instructions:

Utilise the feature of the bolt on feet (fitted in each corner) for fixing into ground, particularly when the product is to remain situated in one permanent location.


OPERATING AUTO YOKE



Instruction: pull handle back (towards main body) to activate auto release.

READY TO RECEIVE

CAUTION: Allow sufficient space around the entire crush body for operating safely, keep fingers and limbs well clear of any moving parts.




Instruction: allow animal to enter the chute

ANIMAL YOKED


CAUTION: Do not enter the area in front of the crush until you are certain the animal is fully secured.

Note: arrow & circle indicate animal travel direction and head position when yoked.



Instruction: pull operating handle forward to activate auto release.


ANIMAL RELEASED



Instruction: push on yoke door to activate the auto yoke, making both doors drop back (angled into the chute).

RETURN TO READY TO RECEIVE


CAUTION: ensure the working area is clear before allowing the next animal to enter the chute .




YOKE APERTURE ADJUSTMENT

CAUTION: ensure the working area is clear before allowing the next animal to enter the

Yoke aperture ranges achievable, 240 - 100 (mm) full height







Instructions: lift the handle attached to the link arm to release from its original locking place. move up or down the 'hair comb like' slanted locking plate as required.


SQUEEZE MECHANISM

CAUTION: Disable squeeze mechanism before letting animal into chute!



Squeeze ranges between 760mm—285mm (unlimited internal width of the chute)






Instruction: push/pull handle to trigger squeeze mechanism.

OPERATING SIDE DOORS


Instruction: lift bolt pivot handle to release (four doors individually operated)



OPERATING FULL LENGTH SIDE RAIL

(removed for entire side access)


Instructions: pull operating handle to release and lift rail away to free, place loose rail away from the crush when not in use.



Option to reverse the side rail (situating the side rail operating handle to lock into the front end main body pillar, adjacent to vet access)


OPERATING VET ACCESS PANELS

Vet access panels feasible of either side, all with capability to remove when access is required.



0


Instruction: remove locking pins.




Instruction: lift (so tube end is free from bushes) and lift away to fully remove the panel.


Note: instructions for removal of the top vet access are the same .

PIVOTING VET ACCESS TOP PANEL



OPERATING RUMP RAIL






Instruction: position rump rail between the ratchet system, behind the animal's backside!

CAUTION: Animal must be inside the chute and yoked before positioning the rump rail.

Slide up to be **tight** behind animal's rear.



RECOMMENDED MAINTENANCE

REGULAR MAINTENANCE CHECKS TO THE PRODUCT WILL GIVE EARLY WARNING OF ANY PROBLEMS, SUCH AS DAMAGE, CORROSION OR MISSING PARTS AND WILL ENSURE YOUR PRODUCT IS RUNNING SMOOTHLY AND SAFELY.

- □ General inspection of the entire product.
- Check operation of all spring-loaded latches, adjustment link arms associated with the side doors and auto yoke. Ensure all moving components are free of debris and are functioning correctly.
- Any damage to steelwork (causing exposure of bare metal) should be made good by rubbing down and applying x1 coat of galvanizing coating and x1 of oil based paint.
- Check all fasteners sets are tight. Ensure all parts are in place before each use.
- □ Remove any obstructions and cut away any vegetation growing on the kit.
- Lubricate all moving parts, in particular latches and locking systems
- □ Simply cleaning down after each use will increase the working life due to the highly corrosive effect of animal's urine and/or faeces.

CONTACT

IAE HEADQUARTERS Mossfield Road Adderley Green Stoke-on-Trent England ST3 5BW

T: +44 (0)1782 339320 F: +44 (0)1782 339321 E: <u>sales@iae.co.uk</u>

Company Reg number: 01448116 Company VAT number: 319060967

TECHNICAL DOCUMENT NO. 1002240-01